



Para 12[®] Extend

Para 12 Extend is a whole blood hematology control which offers complete assay values for the three-part white blood cell differential. It features vials with pierceable caps for autosampling and is available in low, normal and high levels. Para 12 Extend is compatible with Beckman Coulter[®], Erma, Horiba[®] Medical, HTI[®] Diagnostics, Mindray[®], Nihon Kohden[®] and Siemens[®] three-part hematology analyzers.

Features

- + Whole blood three-part differential hematology control
- + Available in three clinically significant levels
- + Plastic vials with pierceable caps for analyzer autosampling
- + Compatible with the most popular three-part differential hematology analyzers
- + 30-day open-vial stability
- + 190-day closed-vial stability

Benefits

- + Helps clinical laboratories meet quality requirements
- + Established assay ranges provide less room for error
- + Long shelf life reduces the number of shipments and associated costs
- + Fewer lot-to-lot correlations save time and resources
- + Access to *STATS*[®], a free interlaboratory quality control program with peer group data comparison and online submission

Ordering information

Level	Description	Catalog#
Low, Normal, High	3 x 2.5 mL	218755
Low, Normal, High	3 X 4.5 mL	218776
Low	6 x 2.5 mL	218993
Normal	6 x 2.5 mL	218758
Normal	6 x 4.5 mL	218995
High	6 x 2.5 mL	218994
Low, Normal, High	18 x 2.5 mL	218750
Low, Normal, High	18 x 4.5 mL	218770